




MIRO · CANVAS™

Quick Start Guide

Installation

- 1 Place on a dry, level surface free from vibration (minimum 30cm W x 45cm D)
- 2 At the back of the Canvas, connect the power supply and either:
 - Screw in the Wi-Fi antenna or Connect an Ethernet cord



- 3 Press the black power button to turn on the Canvas.
The instrument may take a few minutes to power up.
- 4 Use the touchscreen Settings menu to:
 - Connect to a Wi-Fi network
 - Level instrument
 - Run Miro Diagnostics
- 5 You are now ready to start a run! Follow the instructions in the assay-specific protocol to set up reagents and load the cartridge.

Tips & Troubleshooting

- Remain connected to the network whenever possible
Network connectivity is not required to run existing protocols, but is critical for updates and optimal support.
- It is recommended to run Miro Diagnostics weekly
- To open the Canvas lid, secure the base of the instrument with one hand, grasp the front center of the lid with the other hand and lift straight up. Close the lid by pressing down on both front corners.
- Leave the cartridge clamp engaged during sample recovery

Instrument Malfunction	Possible Cause & Solution
On-screen Error: Instrument Not Level	Navigate to Instrument Settings, select "Level Instrument" and follow the instructions on the screen to level the instrument.
Unable to connect the instrument to a Wi-Fi network	Ensure Wi-Fi antenna has been securely screwed into the instrument. Alternatively, you may plug in an Ethernet cord and toggle off Wi-Fi on the Settings menu. Note: the Miro Canvas cannot connect to Wi-Fi networks requiring captive portal login (opening of a page which requires input of user name and/or password).
The Miro Canvas does not successfully start up upon installation (indicated by blue LED indicator bar)	Power may not be flowing to the instrument. Ensure power supply cords are securely attached to the instrument, power brick, and wall outlet, and that the wall outlet is supplying power.

Version 2.0

Miroculus Inc. San Francisco, CA
miroculus.com

For research use only. Not for use in diagnostic procedures.