



A VERSATILE MICROFLUIDIC PLATFORM FOR AUTOMATING COMPLEX BIOLOGICAL AND CHEMICAL PROTOCOLS

Eugenia Carvajal, Eduardo Cervantes, Poornasree Kumar, Winnie Chow, Yu-Hung Chen, Mais J. Jebrail and Foteini Christodoulou


Miroculus, 458 Brannan Street, San Francisco, CA 94107, United States

SUMMARY

Utilizing Miro Technology, Miroculus is developing a platform that allows users to easily build complex protocols for next-generation sequencing (NGS) applications and synthetic chemistry as well as run off-the-shelf cartridges for predefined workflows.


MIRO TECHNOLOGY

Miro Technology utilizes electromechanical forces for manipulating fluids in an automated fashion across a surface of patterned electrodes. This technology integrates key operations to perform a wide range of processes:



NGS: EXOME SEQUENCING APPLICATION

- Miro Technology automates the steps required for Agilent SureSelect XT-HS target enrichment workflow.




- System performance is comparable to benchtop (BT), offering advantages of automation and shorter hands-on processing times. Our initial XT-HS library preparation (Figure 2a) followed by SureSelect Focused Exome capture (Figure 2b) on Miro Technology yielded libraries within 2 - 4.5 % of BT performance across a series of sequencing metrics. With further optimization, performances equal to or higher than BT are expected.




HIGH-THROUGHPUT CHEMICAL SYNTHESIS APPLICATION

- In order to expand the utility of this technology to complex chemical synthesis we designed a feasibility study around palladium-catalyzed carbon-nitrogen bond formation.


Palladium-Catalyzed C-N Reaction:



Pd Precatalysts:



Bases:




Initial nanoscale work: Santanilla et. al. Science 2015, 49.

Feasibility Requirements:

- 4 component mixing
- 8 μ L reaction droplets
- Dispense from stock droplet
- Operate under inert atmosphere
- Quench reactions on chip
- Observe spread in reactivity

- Analysis of the high throughput reactions on our system showed good correlation to previously reported reactivity trends for this catalytic transformation.



Acknowledgements

We thank all of the Team Members within Miroculus and **Discovery High-Throughput Chemistry at GSK** for their support and assistance with projects.